
Traveling wave ion transport for the NSCL cyclotron gas stopper



Maxime Brodeur National Superconducting Cyclotron Laboratory



National Science Foundation Michigan State University

Ion guiding for the cyclotron gas stopper

Cyclotron gas stopper (cycstopper): reversed operation mode cyclotron (S. Schwarz talk)

- Provide thermal RIB to precision experiments (M. Redshaw talk) and for reacceleration (D. Leitner talk).
- Some ions of interest will have $T_{1/2} < 50$ ms and/or produced at low yields.

The efficient and quick transport of the stopped ions is critical.



National Science Foundation Michigan State University


Ion guiding for the cyclotron gas stopper

Cyclotron gas stopper (cycstopper): reversed operation mode cyclotron (S. Schwarz talk)

- Provide thermal RIB to precision experiments (M. Redshaw talk) and for reacceleration (D. Leitner talk).
- Some ions of interest will have $T_{1/2} < 50$ ms and/or produced at low yields.

The efficient and quick transport of the stopped ions is critical.

Performed detailed R&D of transport methods for the cycstopper.



<u>Outline</u>

- 1) Transport in stopping region
 - Ion surfing method experimental results
- 2) Transport in extraction region
 - The ion conveyer experimental results




National Science Foundation Michigan State University

lon guiding in the stopping region

Ion transport choice guided by:

- Space constrain: set by pole piece separation
- Clearance along the path of energetic ions
- Geometry: cylindrical symmetry
- Required axial extraction of the ions





National Science Foundation Michigan State University

lon guiding in the stopping region

Ion transport choice guided by:


- Space constrain: set by pole piece separation
- Clearance along the path of energetic ions
- Geometry: cylindrical symmetry
- Required axial extraction of the ions



Easiest choice:

1) Push ions axially (with ΔV)

2) Transport along repelling surface towards orifice





National Science Foundation Michigan State University

lon guiding in the stopping region

Ion transport choice guided by:

- Space constrain: set by pole piece separation
- Clearance along the path of energetic ions
- Geometry: cylindrical symmetry
- Required axial extraction of the ions





Easiest choice:

1) Push ions axially (with ΔV)

2) Transport along repelling surface towards orifice

RF carpets creates such repelling (Wada-san talk)






Maxime Brodeur EMIS 2012 conference December 4th, 2012



National Science Foundation Michigan State University

Transport methods using RF carpets

Potential gradient








National Science Foundation Michigan State University

Transport methods using RF carpets







National Science Foundation Michigan State University

Ion surfing R&D goals





Require efficient transport of light (A ~ 20) short-lived nuclei (< 50ms) for the cycstopper.

- Test experimentally the limitations of the ion surfing method.
 - (should ultimately achieve faster transport than conventional method)
- Gauge used: the maximum transport speed for efficient transport vs.:
- Pressure (increased pressure gives higher stopping efficiency)
- Push field (larger beam rates, more ionization results larger push towards carpet)
- Ion mass (aim at A ~ 20 and possibly lower)



Ion surfing carpet design






- 1) Maximize resonant *f*_{RF}
- Minimize carpet C by using thin substrate
- 2) Minimize carpet pitch a
- Min. manufacturing gap & electrode width

Also added C in series with carpet to reduce overall C.

Result: $f_{RF} = 6.8$ MHz for a = 375 μ m




Test setup for RF carpet R&D



National Science Foundation Michigan State University

Experimental results vs. simulations



Parameters used:


RF amplitude	75, 51 V	
RF frequency	6.6, 6.8 MHz	
Push field	45, 30 V/cm	
Beam current	1 nA	
Pressure	120, 80 mbar	
Travel distance	30 cm	
Wave speed	60, 22.5 m/s	

Good agreement between experimental results and simulations including hard sphere collisions (IonCool code).



National Science Foundation Michigan State University

Experimental results vs. simulations



Parameters used:

RF amplitude	75, 51 V	
RF frequency	6.6, 6.8 MHz	
Push field	45, 30 V/cm	
Beam current	1 nA	
Pressure	120, 80 mbar	
Travel distance	30 cm	
Wave speed	60, 22.5 m/s	

Good agreement between experimental results and simulations including hard sphere collisions (IonCool code).


Locked mode

 $V_{ion} = V_{wave}$



National Science Foundation Michigan State University

Maximum transport speed





National Science Foundation Michigan State University

Ion transport in extraction region


After leaving stopping region, need quick & efficient transport through magnet yoke.

"Traditional" transport methods:

- RFQ ion guide
- RF funnel

<u>Issues for the cycstopper</u>:

- Magnetic field gradient requires high RF frequencies
- Pressure in region: 2 to 30 mbar limits DC gradient, resulting in slow extraction.





National Science Foundation Michigan State University

Ion transport in extraction region

After leaving stopping region, need quick & efficient transport through magnet yoke.

"Traditional" transport methods:


- RFQ ion guide —
- RF funnel

Issues for the cycstopper:

Magnetic field gradient requires high RF frequencies
Pressure in region: 2 to 30 mbar limits DC gradient, resulting in slow extraction.

Proposed solution: ion conveyer

- Stacked rings
- Traveling wave pushed ions axially
- Weak radial focusing







National Science Foundation Michigan State University Maxime Brodeur EMIS 2012 conference December 4th, 2012

Stacked rings 🔨

Ion conveyer first prototype



- 4 phases traveling wave, 19V max, up to 500 kHz
- Ring structure made on the fly from stockroom material (that's why it looks crooked)
- Rb⁺ ions from surface source




Found good agreement between measured efficiency measurements and simulations including hard sphere collisions (IonCool code)



National Science Foundation Michigan State University

Current & future testings



- Designed ion conveyer
- 8 phases traveling wave
- Prototype have been assembled and ready to test.

- RF carpet with circular stripes
- Improve the RF circuit, $f_{RF} = 8.4$ MHz
- Reached 75 m/s for K at 80 mbar
- Next: transport of Na
- Cycstopper RF carpet will be formed of 6 pie-shaped circular segments





National Science Foundation Michigan State University

Ion transport summary

lon specie	³⁹ K	^{85,87} Rb	²⁴ O (t _{1/2} =65ms)
Max. ion velocity (80 mbar)	~50 m/s	>75 m/s	~10 m/s *
Transport time on carpet *	5 ms	< 4 ms	28 ms
Transport time on conveyer *	4 ms	3 ms	1 ms
Total transport time	9 ms	< 7 ms	29 ms

* Based from simulation results

** Used stopping distribution centroid: 0.27(12) m

- Extraction time lower then 30 ms down to A = 24.
- The transport methods using traveling wave studied will allow both a quick and efficient extraction of radioactive ions.

Special thanks to:

A.E. Gehring, G. Bollen, N. Joshi, S. Schwarz and D.J. Morrissey.



National Science Foundation Michigan State University






National Science Foundation Michigan State University

Ion surfing transport regimes

The travelling wave carries the ions along the carpet in two possible regimes:

- Slipping (if v_{ion} < v_{wave})
- Locked (if v_{ion} = v_{wave})



- With increased wave amp., the ions gets deeper in the trough
- When locked-mode is reached, $v_{ion} = v_{wave}$ remains constant
- Until amplitudes gets too large resulting in ion losses



RF carpet design optimization

For the cycstopper, want efficient transport of light ($A \sim 20$) short-lived nuclei (< 50ms). → Need a carpet with maximized repelling force (since $F_{RF} \propto A$)

$$F_{RF}(\max) = m \cdot 4\pi a f^2 \left(1 - 1.09 \sqrt{\frac{a^2 f}{2 V K_0} \frac{p/p_0}{T/T_0}} \frac{1}{\sin(\pi \gamma 2)} \right)$$

Possible actions: • Increase RF amplitude V (discharge-limited to ~ 75V at 100 mbar)

- Reduce gap/pitch ratio (*F_{RF}* changes weakly with gap/pitch)
- Reduce pitch *a* if
- Increase RF frequency f

Need to reduce carpet capacitance




National Science Foundation Michigan State University

Ion transport using RF & DC gradient

One common transport method makes use of RF carpets *.

1. Made of electrode stripes on which an alternating RF voltage with reversed polarity on adjacent electrode is applied. Results in repealing pseudo-potential.

2. Couples a constant voltage gradient to guide the ions to the extraction orifice.




3. The force from the RF is used to balance a pushing field that drives the ions towards the carpet.


* M. Wada *et al.*, NIM A **204**, 570 (2003)



National Science Foundation Michigan State University




Transport efficiency measurement





National Science Foundation Michigan State University

Transport efficiency measurement



to obtain half-distance x₀




National Science Foundation Michigan State University


Transport velocity measurement

Transport velocity measurement procedure:

Switch down ion source bias (block ions).
 Results in a drop of I_{pads}

3) Travel time: moment drop occurs (obtained using a Boltzmann+sine fit function)







National Science Foundation Michigan State University

Transport velocity measurement

Transport velocity measurement procedure: 34 1) Switch down ion source bias (block ions). 2) Results in a drop of I_{pads} 46' 3) Travel time: moment drop occurs (obtained using a Boltzmann+sine fit function) Rb⁺ or K⁺ scope d = 2.1 cm 0.00 4) Repeat for different distances d *d* = 32.8 cm Picoampmeter output voltage (V) 5) Average ion velocity: $d = \overline{v} \cdot (t - t_{drop})$ -0.02 35 slope: 26.3(6) m/s 30 0.04 distance from pads (cm) 25 0.06 20 15 -0.08 10 5 0.00 0.05 0.10 -0.05 0.15 Time after switching (s) 0 34 42 36 38 40 44 46 48 average total travel time (ms)

S NSCL

National Science Foundation Michigan State University